

Medical Chemical Corp.
19430 Van Ness Ave.
Torrance, CA 90501

Customer Service: Phone (310)787-6800
FAX (310)787-4464

CHEMTREC Emergency Response Telephone Number: (800)424-9300

Note: The CHEMTREC phone number is only for emergencies involving spills, leaks, fire, exposure or accident. Please direct all other inquiries to our customer service phone number.

Section I - Product Identification

A mixture of terphenyl, hydrogenated terphenyl, polybutenes and natural hydrocarbons.

Section II - Hazards Identification

Overview: May be harmful if swallowed. May be irritating to skin eyes and respiratory tract.

Safety Ratings


Health: Normal Flammability: Slight Reactivity: None Contact: Slight

Recommended safety equipment: safety goggles, lab coat and proper gloves

Storage: General storage

NFPA Ratings

Health = 0 Flammability = 1 Reactivity = 0



Potential Health Effects

The toxicology of this compound have not been completely examined. It is presumed that the toxicity of this item is similar to other hydrocarbons.

Inhalation: May be irritating

Ingestion: While the toxicity of this compound is low, large doses may cause nausea, vomiting, diarrhea, etc.

Skin contact: Not normally a problem

Eye contact: May be irritating

Chronic Exposure: Unknown

Aggravation of preexisting conditions: Unknown

Section III - Composition/Information on Components

Ingredients	CAS#	OSHA Pel	ACGIH TLV	Other Limits	%
Hydrogenated terphenyl	61788-32-7	5 mg/m3	5 mg/m3		trade secret
Terphenyl	92-94-4	5 mg/m3	5 mg/m3		trade secret
Polybutenes	—	5 mg/m3	5 mg/m3		trade secret
Natural hydrocarbons	—	5 mg/m3	5 mg/m3		trade secret

Section IV - First Aid Measures

Inhalation: Unlikely to be a problem under normal circumstances. Remove from source of exposure and get medical attention for any breathing difficulty.

Ingestion: Do not induce vomiting. Never give anything by mouth to an unconscious person.

Skin Contact: Wash affected area with soap and water. Get medical advice if irritation develops.

Eye Contact: Rinse thoroughly with running water. Get medical advice if irritation develops.

Section V - Fire Fighting Measures

Flash point: 163 °C (325 °F) TOC

Flammable Limits: Not applicable.

Fire: Not normally a fire Hazard.

Explosion: Not Normally an explosion hazards.

Fire Extinguishing Media: Any means suitable for surrounding fire.

Special information: Pyrolysis will release corrosive oxides.

Section VI - Accidental Release Measures

Absorb with a suitable absorbent (such as paper towels) and store in a suitable container for disposal.

Section VII - Handling and Storage

Store in a closed container, protected from freezing.

Section VIII - Exposure Control/Personal Protection

Airborne Exposure Limits: See section III.

Ventilation System: Usually not required. When required, Refer to the ACGIH document, "Industrial Ventilation, a Manual of Recommended Practices" for details about ventilation.

Personal Respirator: Usually not required. In case of emergency, or when exposure levels are unknown, use a positive pressure, full face piece, air supplied respirator.

Skin protection: Protective gloves are not required but recommended as part of good laboratory practice.

Eye Protection: Laboratory safety goggles or similar products are not required but recommended as part of good laboratory practice.

Section IX - Physical and Chemical Properties

Boiling Point: ~340 °C (644 °F)

Vapor pressure (mm Hg): < 0.1 @ 20°C

Vapor Density (air = 1): ~1

Appearance and Odor: A clear colorless liquid.

Density: 0.9 g/ml

Evaporation Rate (water = 1): 0

Solubility: Insoluble

Section X - Stability and Reactivity

Stability: Freezes at very low temperature.

Hazardous Decomposition Products: Nothing unusual.

Hazardous polymerization: Will not occur.

Incompatibilities: Nothing unusual.

Conditions to avoid: Excessive cold/heat and light.

Section XI - Toxicological Information

None relating to normal exposure.

Cancer lists

<u>Ingredient</u>	<u>Known Carcinogenicity?</u>	<u>NTP?</u>	<u>Anticipated?</u>	<u>IARC Category</u>
Immersion Oil	No	No	No	None

Section XII - Ecological Information

Environmental Fate: Unknown.

Environmental Toxicity: Unknown, but not expected to be a problem.

Section XIII - Disposal Considerations

Usually not restricted but Local governments can restrict the amounts of hydrocarbons that may be flushed down drain. Insure compliance with all government regulations.

Section XIV - Transportation Information

Not regulated.

Section XV - Regulatory Information**Chemical Inventory Status**

<u>Ingredient</u>	<u>TSCA</u>	<u>EC</u>
Immersion Oil	Yes	Yes

Federal, State and International Regulations

<u>Ingredient</u>	<u>SARA 302</u>		<u>SARA 313</u>		<u>RCRA</u>	<u>TSCA</u>	<u>Ca. Prop 65</u>
	<u>RQ</u>	<u>TPQ</u>	<u>List</u>	<u>Category</u>	<u>261.33</u>	<u>8(D)</u>	
Immersion Oil	No	No	No	No	No	No	No

Chemical Weapons Convention: No TSCA 12(b): No CDTA: No

SARA 311/312: Acute: No, Chronic: No, Fire: Yes

Section XVI - Other Information

This information is believed to be correct but is not warranted as such, nor does it purport to be all inclusive.

Revision Date: Jun. 25, 2019