

Creation date: 04-21-14

Revision date: 03-01-21

Revision number: 1.0

Section I - Product Identification

Product: An aqueous solution of Potassium hydrogen phthalate.

Intended Uses: A reducing reagent for general laboratory use.

Uses advised against: Not for food use. Intended for laboratory use only.


Manufacturer Identification

Medical Chemical Corp.
19430 Van Ness Ave.
Torrance, CA 90501

Customer Service: Phone (310)787-6800
Email: Christinaavena@med-chem.com
FAX (310)787-4464

Emergency Telephone Number

CHEMTREC Emergency Response Telephone Number: (800)424-9300. Note: The CHEMTREC phone number is only for emergencies involving spills, leaks, fire, exposure or accident. Please direct all other inquiries to our customer service phone number.



Section II - Hazard Identification

This item is not considered hazardous by 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200).

Acute toxicity (Oral): Category 5 (H303). May be harmful if swallowed.

Acute toxicity (Dermal): Category 5 (H313). May be harmful in contact with skin.

Eye damage/eye irritation: Category 2B (H320). May cause eye irritation.

The manufacturer is unaware of any chemical hazard associated with this product. The product may cause skin or eye irritation and may be harmful if swallowed.

Precautionary Statements

P102 Keep out of reach of children.
P103 Read label before use.
P270 Do not eat, drink or smoke when using this product.

Safety Ratings

Health: Slight Flammability: None *Reactivity:* Stable *Contact:* None
Recommended safety equipment: Not normally required.
Storage: General storage

NFPA Ratings

Health = 1 Flammability = 0 Reactivity = 0

Potential Health Effects

The manufacturer is not aware of any serious hazards associated with this product.

Inhalation: May be Irritating to the respiratory tract.

Ingestion: While the toxicity of this compound is low, large doses may cause nausea, vomiting, diarrhea, etc

Skin contact: Can cause skin irritation or aggravation of existing dermatitis.

Eye contact: Can cause eye irritation.

Chronic Exposure: Not likely to be an issue.

Aggravation of preexisting conditions: Unknown.

Section III - Composition/Information on Components

| Ingredient | CAS# | EC/List no. | % w/w |
|------------------------------|----------|-------------|-------|
| Potassium hydrogen phthalate | 877-24-7 | 212-889-4 | 1% |

Also contains trace amounts of a preservative.

Section IV - First Aid Measures

Inhalation: Unlikely to be an issue.

Ingestion: Unlikely to be an issue. Get medical advise for exceptional cases.

Skin Contact: Unlikely to be an issue. Wash affected area with soap and water. Get medical advice if irritation develops.

Eye Contact: Rinse thoroughly with running water. Get medical advice if irritation develops.

Section V - Fire Fighting Measures

Flash point: Not applicable.

Flammable Limits: Not applicable.

Fire: Not normally a fire hazard.

Explosion: Not normally an explosion hazard.

Fire Extinguishing Media: Any means suitable for surrounding fire.

Special information: Nothing exceptional.

Section VI - Accidental Release Measures

Absorb with a suitable absorbent and store in a suitable container for disposal.

Section VII - Handling and Storage

P403+P233+P102; Store in a well-ventilated place. Keep container tightly closed. Keep out of reach of children.

Section VIII - Exposure Control/Personal Protection

| Component | ACGIH TLV | OSHA PEL | NIOSH IDLH |
|------------------------------|---------------|---------------|---------------|
| Potassium hydrogen phthalate | Not regulated | Not regulated | Not regulated |

Legend

ACGIH: American Conference of Governmental Industrial Hygienists.

OSHA: Occupational Safety and Health Administration.

NIOSH: National Institute for Occupational Safety and Health.

IDLH: Immediately dangerous to life or health.

Ventilation System: Usually not required. When required, Refer to the ACGIH document, "Industrial Ventilation, a Manual of Recommended Practices" for details about ventilation.

Personal Respirator: Usually not required. In case of emergency, or when exposure levels are unknown, use a positive pressure, full face piece, air supplied respirator.

Skin protection: Protective gloves are not required but may be used as part of good laboratory practice.

Eye Protection: Laboratory safety goggles or similar products are not required but recommended as part of good laboratory practices

Section IX - Physical and Chemical Properties

Boiling Point: 100 °C

Vapor pressure (mm Hg): 18 @ 20 °C

Vapor Density (air = 1): 0.6

Volatile organic carbon (VOC): Not applicable

Density: 1.00 g/ml

Evaporation Rate (water = 1): 1

Solubility: Infinitely miscible with water

Appearance and Odor: A clear, colorless and odorless reagent.

Section X - Stability and Reactivity

Stability: Stable under normal conditions.

Hazardous Decomposition Products: No known hazardous decomposition products.

Hazardous polymerization: Will not occur.

Incompatibilities: Nothing unusual.

Conditions to avoid: Nothing unusual.

Section XI - Toxicological Information

| Component | LD50 Oral | LD50 Dermal | LC50 Inhalation |
|------------------------------|-------------------|-------------|-----------------|
| Potassium hydrogen phthalate | >3200 mg/kg (Rat) | No data | No data |

The product is non-toxic by ingestion or inhalation.

Cancer lists

| <i>Ingredient</i> | <i>Known Carcinogenicity?</i> | <i>NTP?</i> | <i>Anticipated?</i> | <i>IARC Category</i> |
|------------------------------|-------------------------------|-------------|---------------------|----------------------|
| Potassium hydrogen phthalate | No | No | No | None |

Section XII - Ecological Information

This product is not expected to bioaccumulate. The material is removed from the air by dry and liquid adsorption. The product biodegrades quickly and does not contain endocrine disruptors, persistent organic pollutants or ozone depleters.

Soil Mobility: Unknown

Environmental Fate: Biodegradable

Environmental Toxicity: The product has low environmental toxicity.

| Component | Freshwater Fish | Water Flea | Freshwater algae |
|------------------------------|-----------------|------------|------------------|
| Potassium hydrogen phthalate | No data | No data | No data |

Section XIII - Disposal Considerations

Not normally restricted but ultimate authority resides with local government agencies. Insure compliance with all government regulations.

Section XIV - Transportation Information

Not regulated by DOT or IATA.

Section XV - Regulatory Information**Chemical Inventory Status**

| <i>Ingredient</i> | <i>TSCA</i> | <i>EC</i> |
|------------------------------|-------------|-----------|
| Potassium hydrogen phthalate | Yes | Yes |

Federal, State and International Regulations

| <i>Ingredient</i> | <i>SARA 302</i> | | <i>SARA 313</i> | | <i>RCRA</i> | <i>TSCA</i> | |
|------------------------------|-----------------|------------|-----------------|-----------------|---------------|-------------|--------------------|
| | <i>RQ</i> | <i>TPQ</i> | <i>List</i> | <i>Category</i> | <i>261.33</i> | <i>8(D)</i> | <i>Ca. Prop 65</i> |
| Potassium hydrogen phthalate | No | No | No | No | No | No | No |

Chemical Weapons Convention: No TSCA 12(b): No CDTA: No

SARA 311/312: *Acute*: No, *Chronic*: No, *Fire*: No, *Pressure*: No

Section XVI - Other Information

This information is believed to be correct at the time of publication but is not guaranteed as such, nor does it purport to be all inclusive. Medical Chemical Corp. assumes no liability for the accuracy or completeness of the information. The user assumes all responsibility for compliance with federal, state and local laws.

Document Number: S108

Revision number: 1.0

Revision Date: Mar. 1, 2021