

Creation date: 04-01-15

Revision date: 04-01-21

Revision number: 1.0

Section I - Product Identification

Product: ACS grade deionized water. Meets the specifications for ASTM type I water.

Intended Uses: A general purpose laboratory reagent.

Uses advised against: Not for food use. Intended for laboratory use only.


Manufacturer Identification

Medical Chemical Corp.
19430 Van Ness Ave.
Torrance, CA 90501

Customer Service: Phone (310)787-6800
Email: Christinaavena@med-chem.com
FAX (310)787-4464

Emergency Telephone Number

CHEMTREC Emergency Response Telephone Number: (800)424-9300. Note: The CHEMTREC phone number is only for emergencies involving spills, leaks, fire, exposure or accident. Please direct all other inquiries to our customer service phone number.



Section II - Hazard Identification

This item is not considered hazardous by 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200).

Acute toxicity (Oral): Category 5 (H303). May be harmful if swallowed.

Acute toxicity (Dermal): Category 5 (H313). May be harmful in contact with skin.

Eye damage/eye irritation: Category 2B (H320). May cause eye irritation.

The manufacturer is unaware of any chemical hazard associated with this product. The product may be harmful if swallowed in excess.

Precautionary Statements

- P102 Keep out of reach of children.
- P103 Read label before use.
- P270 Do not eat, drink or smoke when using this product.

Safety Ratings

Health: None *Flammability:* None *Reactivity:* None *Contact:* None
Recommended safety equipment: Not normally required.

NFPA Ratings

Health = 0 Flammability = 0 Reactivity = 0

Potential Health Effects

The manufacturer is not aware of any serious hazards associated with this product..

Inhalation: May be Irritating to the respiratory tract.

Ingestion: While the toxicity of this compound is low, large doses may cause nausea, vomiting, diarrhea, etc

Skin contact: Can cause skin irritation or aggravation of existing dermatitis.

Eye contact: Can cause eye irritation.

Chronic Exposure: Not likely to be an issue.

Aggravation of preexisting conditions: Unknown.

Section III - Composition/Information on Components

Ingredient	CAS#	EC/List No.	%
Water	7732-18-5	231-791-2	100.0% w/w

Section IV - First Aid Measures

Inhalation: Unlikely to be an issue.

Ingestion: Unlikely to be an issue. Hey, it's just water!

Skin Contact: Unlikely to be an issue. Get medical advice if irritation develops.

Eye Contact: Rinse thoroughly with running water. Get medical advice if irritation develops.

Accidental immersion: Summon a lifeguard.

Section V - Fire Fighting Measures

Flash point: Not applicable.

Flammable Limits: Not applicable.

Fire: Extremely unlikely to be a fire hazard.

Explosion: Not normally an explosion hazard.

Fire Extinguishing Media: Any means suitable for the surrounding fire.

Special information: Nothing exceptional.

Section VI - Accidental Release Measures

Absorb with a suitable absorbent and store in a suitable container for disposal.

Section VII - Handling and Storage

Store in a closed container at room temperature.

Section VIII - Exposure Control/Personal Protection

Component	ACGIH TLV	OSHA PEL	NIOSH IDLH
Water	Not regulated	Not regulated	Not regulated

Legend

ACGIH: American Conference of Governmental Industrial Hygienists.

OSHA: Occupational Safety and Health Administration.

NIOSH: National Institute for Occupational Safety and Health.

IDLH: Immediately dangerous to life or health.

Ventilation System: Usually not required. When required, Refer to the ACGIH document, "Industrial Ventilation, a Manual of Recommended Practices" for details about ventilation.

Personal Respirator: Usually not required. In case of emergency, or when exposure levels are unknown, use a positive pressure, full face piece, air supplied respirator.

Skin protection: Protective gloves are not required but may be used as part of good laboratory practice.

Eye Protection: Laboratory safety goggles or similar products are not required but recommended as part of good laboratory practices

Section IX - Physical and Chemical Properties

Boiling Point: 100 °C

Vapor pressure (mm Hg): 18 @ 20 °C

Vapor Density (air = 1): 0.6

Appearance and Odor: A clear, colorless and odorless Solution.

Density: About 0.99766 g/ml @ 22.5 °C

Evaporation Rate (water = 1): 1

Solubility: Infinitely miscible with water

Section X - Stability and Reactivity

Stability: Stable under normal conditions.

Hazardous Decomposition Products: No known hazardous decomposition products.

Hazardous polymerization: Will not occur.

Incompatibilities: Anything that reacts with water.

Conditions to avoid: Nothing unusual.

Section XI - Toxicological Information

Component	LD50 Oral	LD50 Dermal	LC50 Inhalation
Water	No data	No data	No data

The product is non-toxic by ingestion or inhalation.

Cancer lists

<i>Ingredient</i>	<i>Known Carcinogenicity?</i>	<i>NTP?</i>	<i>Anticipated?</i>	<i>IARC Category</i>
Water	No	No	Not hardly	None

Section XII - Ecological Information

This product is not expected to bioaccumulate. The material is removed from the air by dry and liquid adsorption. The product does not contain endocrine disrupters, persistent organic pollutants or ozone depleters.

Soil Mobility: Extremely mobile in soil.

Environmental Fate: Biodegradable

Environmental Toxicity: The product has no environmental toxicity except in very large quantities.

Component	Freshwater Fish	Water Flea	Freshwater algae
Water	No data	No data	No data

Section XIII - Disposal Considerations

At the present time there is no neutralizer for this product. Not restricted but local governments can restrict the amounts of anything that can be flushed down drain. Dispose of contents and container in accordance with all government regulations.

Section XIV - Transportation Information

Not regulated by DOT or IATA.

Section XV - Regulatory Information

Chemical Inventory Status

<i>Ingredient</i>	<i>TSCA</i>	<i>EC</i>
Water	Yes	Yes

Federal and State Regulations

<i>Ingredient</i>	<i>SARA 302 RQ</i>	<i>SARA 302 TPQ</i>	<i>SARA 313 List</i>	<i>SARA 313 Category</i>	<i>RCRA 261.33</i>	<i>TSCA 8(D)</i>	<i>Ca. Prop 65</i>
Water	No	No	No	No	No	No	No

Chemical Weapons Convention: No TSCA 12(b): No CDTA: No

SARA 311/312: *Acute:* No, *Chronic:* No, *Fire:* No, *Pressure:* No

Section XVI - Other Information

This information is believed to be correct at the time of publication but is not guaranteed as such, nor does it purport to be all inclusive. Medical Chemical Corp. assumes no liability for the accuracy or completeness of the information. The user assumes all responsibility for compliance with federal, state and local laws.

Document Number: S049

Revision number: 1.0

Revision Date: Apr. 1, 2021