

Creation date: 05-01-15

Revision date: 03-30-21

Revision number: 1.2

Section I - Product Identification

Product: An aqueous solution of povidone iodine.

Intended Uses: An in-vitro diagnostic reagent used for gram staining. Also used as a general laboratory reagent.

Uses advised against: Not intended to be used as an antiseptic.

Manufacturer Identification

Medical Chemical Corp.
19430 Van Ness Ave.
Torrance, CA 90501

Customer Service: Phone (310)787-6800
Email: Christinaavena@med-chem.com
FAX (310)787-4464

Emergency Telephone Number

CHEMTREC Emergency Response Telephone Number: (800)424-9300. Note: The CHEMTREC phone number is only for emergencies involving spills, leaks, fire, exposure or accident. Please direct all other inquiries to our customer service phone number.



Section II - Hazard Identification

This item is not considered hazardous by 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200).

Eye damage/eye irritation: Category 2 (H319). May cause serious eye irritation.

Acute toxicity (Oral): Category 4 (H302). Harmful if swallowed.

Acute toxicity (Dermal): Category 4 (H312). Harmful in contact with skin

Signal word: Warning

Hazard Statements: May cause skin and eye irritation. Harmful if swallowed.

Precautionary Statements

P102 Keep out of the reach of children.

Safety Ratings

Health: Slightly Hazardous *Flammability:* None *Reactivity:* Unstable if heated *Contact:* Slight
Recommended safety equipment: safety goggles, lab coat and proper gloves are recommended as part of good laboratory practice.

NFPA Ratings

Health = 2 Flammability = 0 Reactivity = 1

Potential Health Effects

The toxicology of this compound has not been completely examined. It is presumed that the toxicity of this item is similar to other dilute Iodine solutions.

Inhalation: May be Irritating to respiratory tract.

Ingestion: While the toxicity of this compound is low, large doses may cause nausea, vomiting, diarrhea, etc

Skin contact: Not normally a problem. May cause skin irritation or aggravation of existing dermatitis.

Eye contact: May cause serious eye irritation.

Chronic Exposure: Unknown.

Aggravation of preexisting conditions: Unknown.

Section III - Composition/Information on Components

Ingredient	CAS #	EC/List No.	%
Povidone Iodine	25655-41-8	924-458-4	5%-15% w/w

Section IV - First Aid Measures

Inhalation: Remove from source of exposure and get medical attention for any breathing difficulty.

Ingestion: Do not induce vomiting. Aspiration into the lungs may produce death. Get immediate medical attention even if symptoms improve.

Skin Contact: Unlikely to be an issue. Wash affected area with soap and water. Get medical advice if irritation develops.

Eye Contact: Rinse thoroughly with running water. Get medical advice if irritation develops.

Section V - Fire Fighting Measures

Flash point: Not applicable.

Flammable Limits: Not applicable.

Fire: Not normally a fire hazard.

Explosion: Not normally an explosion hazard.

Fire Extinguishing Media: Any means suitable for surrounding fire.

Special information: Pyrolysis will release corrosive Iodine vapor.

Section VI - Accidental Release Measures

Absorb with a suitable absorbent and store in a suitable container for disposal.

Section VII - Handling and Storage

P403+P233+P102; Store in a well-ventilated place. Keep container tightly closed. Keep out of reach of children. Iodine is light sensitive.

Section VIII - Exposure Control/Personal Protection

Component	ACGIH TLV	OSHA PEL	NIOSH IDLH
Povidone Iodine	Not Established	Not Established	Not Established

Legend

ACGIH: American Conference of Governmental Industrial Hygienists.

OSHA: Occupational Safety and Health Administration.

NIOSH: National Institute for Occupational Safety and Health.

IDLH: Immediately dangerous to life or health.

Ventilation System: Usually not required. When required, Refer to the ACGIH document, "Industrial Ventilation, a Manual of Recommended Practices" for details about ventilation.

Personal Respirator: Usually not required. In case of emergency, or when exposure levels are unknown, use a positive pressure, full face piece, air supplied respirator.

Skin protection: Protective gloves are not required but may be used as part of good laboratory practice.

Eye Protection: Laboratory safety goggles or similar products are not required but recommended as part of good laboratory practices

Section IX - Physical and Chemical Properties

Boiling Point: 97 °C

Vapor Pressure (mm Hg): 18 @ 20 °C

Vapor Density (air = 1): 0.6

Appearance and Odor: An opaque purple solution with an acrid odor.

Density: About 1.025 g/ml

Evaporation Rate (water = 1): 1

Solubility: Infinitely miscible with water

Section X - Stability and Reactivity

Stability: Stable under normal conditions.

Hazardous Decomposition Products: Decomposes to form elemental Iodine.

Hazardous polymerization: Will not occur.

Incompatibilities: Nothing unusual.

Conditions to avoid: Heat and light.

Section XI - Toxicological Information

Component	LD50 Oral	LD50 Dermal	LC50 Inhalation
Povidone Iodine	8100 mg/kg (Rat)	Not Available	Not Available

Cancer lists

<i>Ingredient</i>	<i>Known Carcinogenicity?</i>	<i>NTP?</i>	<i>Anticipated?</i>	<i>IARC Category</i>
Povidone Iodine	No	No	No	None

Section XII - Ecological Information

This product is not expected to bioaccumulate. The material is removed from the air by dry and liquid adsorption. The product biodegrades quickly and does not contain endocrine disruptors, persistent organic pollutants or ozone depleters.

Soil Mobility: Unknown

Environmental Fate: Biodegradable

Environmental Toxicity: The product has low environmental toxicity

Component	Freshwater Fish	Water Flea	Freshwater algae
Povidone Iodine	Not Available	Not Available	Not Available

Section XIII - Disposal Considerations

Not normally restricted but ultimate authority resides with local government agencies. Insure compliance with all government regulations

Section XIV - Transportation Information

Not regulated by DOT or IATA.

Section XV - Regulatory Information

Chemical Inventory Status

<i>Ingredient</i>	<i>TSCA</i>	<i>EC</i>
Povidone Iodine	Yes	Yes

Federal and State Regulations

<i>Ingredient</i>	<i>SARA 302</i>		<i>SARA 313</i>		<i>RCRA</i>	<i>TSCA</i>	<i>Ca. Prop 65</i>
	<i>RQ</i>	<i>TPQ</i>	<i>List</i>	<i>Category</i>	<i>261.33</i>	<i>8(D)</i>	
Povidone Iodine	No	No	No	No	No	No	No

Chemical Weapons Convention: No TSCA 12(b): No CDTA: No

SARA 311/312: *Acute:* Yes, *Chronic:* No, *Fire:* No, *Pressure:* No

Section XVI - Other Information

This information is believed to be correct at the time of publication but is not guaranteed as such, nor does it purport to be all inclusive. Medical Chemical Corp. assumes no liability for the accuracy or completeness of the information. The user assumes all responsibility for compliance with federal, state and local laws.

Document Number: S151

Revision number: 1.2

Revision Date: Mar. 30, 2021