

Creation date: 04-15-2021

Revision date: 09-13-2024

Revision number: 1.05

Section I - Product Identification

Product: A dilute aqueous solution of Periodic acid.

Intended Uses: An in-vitro diagnostic reagent used with Schiff's stain. Also used as a general laboratory reagent.

Uses advised against: For laboratory use only.

Country of origin: United States.

Manufacturer Identification

Medical Chemical Corp.
19430 Van Ness Ave.
Torrance, CA 90501

Customer Service: Phone (310)787-6800
Email. Customerservice@med-chem.com
FAX (310)787-4464

Emergency Telephone Number

CHEMTREC Emergency Response Telephone Number: (800)424-9300. Note: The CHEMTREC phone number is only for emergencies involving spills, leaks, fire, exposure or accident. Please direct all other inquiries to our customer service phone number.



Section II - Hazard Identification

This item is not considered hazardous by 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200).

Eye damage/eye irritation: Category 2 (H319). May cause serious eye irritation.

Acute toxicity (Oral): Category 4 (H302). Harmful if swallowed.

Acute toxicity (Dermal): Category 4 (H312). Harmful in contact with skin

Corrosive to Metals: Category 1 (H290), May be corrosive to metals.

Signal word: Warning

- P102 Keep out of reach of children.
P220 Keep away from clothing and other combustible materials.
P260 Do not breathe mist, vapors or spray.
P264 Wash thoroughly after handling.
P270 Do not eat, drink or smoke when using this product.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P363 Wash contaminated clothing before reuse.

Safety Ratings

Health: Slight Flammability: None Reactivity: Stable Contact: Slight
Recommended safety equipment: Not normally required.

NFPA Ratings

Health = 1. Can cause significant irritation.
Flammability = 0. Will not burn
Reactivity = 0. Normally stable, even under fire conditions.

Potential Health Effects

The toxicology of this compound have not been completely examined. It is presumed that the toxicity of this item is similar to other solutions of dilute acids.

Inhalation: May be Irritating to respiratory tract.

Ingestion: While the toxicity of this compound is low, large doses will cause nausea, vomiting, diarrhea, etc

Skin contact: May cause skin damage. May cause skin irritation or aggravation of existing dermatitis.

Eye contact: A serious eye irritant.

Chronic Exposure: Unknown.

Aggravation of preexisting conditions: Unknown.

Section III - Composition/Information on Components

Ingredient	CAS#	EC/List-No.	% w/w
Periodic acid	13444-71-8	236-585-6	1%

Section IV - First Aid Measures

Inhalation: Unlikely to be an issue. Remove from source of exposure and get medical attention for any breathing difficulty.

Ingestion: Do not induce vomiting. Contact a poison control center immediately

Skin Contact: Wash affected area with soap and water. Get medical advice if irritation develops.

Eye Contact: Rinse thoroughly with running water. Get medical advice.

Section V - Fire Fighting Measures

Fire: Not normally a fire hazard.

Explosion: Not normally an explosion hazard.

Fire Extinguishing Media: Any means suitable for surrounding fire.

Special information: Not applicable.

Section VI - Accidental Release Measures

Absorb with a suitable absorbent such as paper towels and store in a suitable container for disposal.

Section VII - Handling and Storage

P403+P233+P102; Store in a well-ventilated place. Keep container tightly closed. Keep out of reach of children.

Section VIII - Exposure Control/Personal Protection

Component	ACGIH TLV	OSHA PEL	NIOSH IDLH
Periodic acid	Not listed	Not listed	Not listed

Legend

ACGIH: American Conference of Governmental Industrial Hygienists.

OSHA: Occupational Safety and Health Administration.

NIOSH: National Institute for Occupational Safety and Health.

IDLH: Immediately dangerous to life or health.

Ventilation System: Usually not required. When required, Refer to the ACGIH document, "Industrial Ventilation, a Manual of Recommended Practices" for details about ventilation.

Personal Respirator: Usually not required. In case of emergency, or when exposure levels are unknown, use a positive pressure, full face piece, air supplied respirator.

Skin protection: Protective gloves are not required but may be used as part of good laboratory practice.

Eye Protection: Laboratory safety goggles or similar products are not required but recommended as part of good laboratory practices

Section IX - Physical and Chemical Properties

Appearance and Odor: A clear, colorless and odorless solution.

Flash point: Not applicable.

Auto ignition temperature: Not applicable

Boiling point range: 99 - 101 °C

Decomposition temperature: No data

Melting point: 0.0 °C

pH: No data

Vapor pressure (mm Hg): 18 @ 20 °C

Viscosity: ~ 0.01 poise

Flammable Limits: Not applicable.

Boiling Point: 100 °C (212 °F)

Density: 1.00 g/ml @ 25 °C

Evaporation Rate (Water = 1): 1

Partition coefficient (n-Octanol/Water): No data

Solubility: Infinitely miscible with water

Vapor density (air = 1): 0.6

Volatile organic carbon (VOC): Not applicable

Section X - Stability and Reactivity

Stability: Stable under normal conditions.

Hazardous Decomposition Products: No known hazardous decomposition products.

Hazardous polymerization: Will not occur.

Incompatibilities: Dimethylsulfoxide, strong bases and strong acids.

Conditions to avoid: Nothing unusual.

Section XI - Toxicological Information

Component	LD50 Oral	LD50 Dermal	LC50 Inhalation
Periodic acid	No data	No data	No data

Cancer lists

Ingredient *Known Carcinogenicity?* *NTP?* *Anticipated?* *IARC Category*

Periodic acid No No No No

Aspiration hazard: No data.

Developmental effects: No data

Irritation: No data

Mutagenic effects: No data

Sensitization: No data

Section XII - Ecological Information

This product is not expected to bioaccumulate. The product biodegrades quickly and does not contain endocrine disrupters, persistent organic pollutants or ozone depleters.

Soil Mobility: Unknown

Environmental Fate: Biodegradable

Environmental Toxicity: The product has low environmental toxicity.

Component	Freshwater Fish	Water Flea	Freshwater algae
Periodic acid	No data	No data	No data

Section XIII - Disposal Considerations

Not normally restricted but ultimate authority resides with local government agencies. In jurisdictions that restrict drain disposal of dilute acids the usual rule is that the sewage exiting the building can't have a pH below 4. Insure compliance with all government regulations. Dilute acids are corrosive to metal plumbing

Section XIV - Transportation Information

Not regulated by DOT or IATA.

Section XV - Regulatory Information

Chemical Inventory Status

Ingredient *TSCA* *EC*

Periodic acid Yes Yes

Federal and State Regulations

Ingredient *SARA 302* *SARA 313* *RCRA* *TSCA*
RQ *TPQ* *List* *Category* *261.33* *8(D)* *Ca. Prop 65*

Periodic acid No No No No No No

Chemical Weapons Convention: No TSCA 12(b): No CDTA: No

SARA 311/312

Acute: Yes

Chronic: No

Fire: No

Pressure: No

Section XVI - Other Information

This information is believed to be correct at the time of publication but is not guaranteed as such, nor does it purport to be all inclusive. Medical Chemical Corp. assumes no liability for the accuracy or completeness of the information. The user assumes all responsibility for compliance with federal, state and local laws.

Document Number: S133

Revision number: 1.05

Revision Date: Sep. 13, 2024